

Clinicare

Impression Disinfectant

Instrument Grade - Low Level



Trigger Spray

Reordering Code: DL2250

- **PREVENTS CROSS-CONTAMINATION**

protect staff and dental technicians from exposure to saliva and blood-borne pathogens.

- **MEETS INFECTION CONTROL GUIDELINES**

complies with ADA guidelines for impression disinfection.

- **SAFE ON IMPRESSION MATERIALS**

non-corrosive formula, no phenol or formaldehyde.

- **QUICK AND EFFECTIVE**

broad-spectrum disinfection in 1 minute.

- **SUPPORTS BEST PRACTICE**

essential step in maintaining a clean, compliant, and professional dental workflow.

Effective Against (contact/kill time) TESTED UNDER DIRTY CONDITIONS

Enterococcus hirae	1 minute
Pseudomonas aeruginosa	1 minute
Staphylococcus aureus	1 minute
Methicillin Resistant Staphylococcus aureus (MRSA)	1 minute
Candida albicans	1 minute
Adenovirus	1 minute
Human herpesvirus 1	1 minute
Norovirus	1 minute



Manufactured by:
Dentalife Australia Pty Ltd
Factory 9/505 Maroondah Highway
Ringwood, Victoria 3134 Australia
Developed and manufactured in Australia
www.dentalife.com.au



BACTERICIDAL YEASTICIDAL