

Clinicare

Medical Enzyme Cleaner



Reordering Code: DL2201

EFFECTIVE FORMULATION WITH MULTIPLE ENZYMES

allows for rapid break down and removes complex organic soils.

SUPPORTS INFECTION CONTROL

complies with ADA guidelines for reprocessing of reusable medical devices when combined with physical cleaning, solutions provides an effective method for cleaning.

LOW FOAMING

is ideal for manual and ultrasonic cleaning

NEUTRAL PH

safe for users and compatible with instruments



Manufactured by:
Dentalife Australia Pty Ltd
Factory 9/505 Maroondah Highway
Ringwood, Victoria 3134 Australia
Developed and manufactured in Australia
www.dentalife.com.au

