






**MK-DENT**






**PROMOTION**

**ADAM**  **DENTAL**  
supplies

## PRIME LINE TURBINES

-  AeroAcoustic blade
-  Reverse flow protection
-  Refined ceramic bearings

-  Glass rod fiber optics
-  Automatic pressure regulation<sup>3</sup>
-  Scratch resistance  
Iconic Prime body

25  
W

Powerful & Precise  
Power Head - 25W

55  
dB

Incredibly Quiet  
Only 55 dB noise



Available for  
KaVo<sup>®</sup>, NSK<sup>®</sup>, W&H<sup>®</sup> and Sirona<sup>®</sup> Couplings



Coupling  
Connection

Kavo  
NSK  
W&H  
Sirona

Power Head  
25 Watt

HP21KL  
HP21NL  
HP21WL  
HP21SL

Small Head  
18 Watt

HP22KL  
HP22NL  
HP22WL  
- N/A -



**BUY 2  
GET 1  
FREE**  
Mix & Match

**\$1,683**





**2 YEARS  
UNCONDITIONAL  
WARRANTY\***

\*Excluding physical damage on  
Prime Line handpieces

## ECO LINE TURBINES



**1 Year FairCare Warranty<sup>1</sup>**  
Unconditional warranty  
except physical damage

-  AeroPerformance blade
-  Refined ceramic bearings
-  Glass rod fiber optics
-  Smooth Eco body

Now with  
24W<sup>2</sup>



Available for  
KaVo<sup>®</sup>, NSK<sup>®</sup>, W&H<sup>®</sup>  
and Sirona<sup>®</sup> Couplings



**BUY 2  
\$890  
each**  
Mix & Match

**2 x Ecoline F/Optic**



**BUY 3+  
\$792  
each**  
Mix & Match

**3 x Ecoline F/Optic**

Coupling  
Connection

KaVo  
NSK  
W&H  
Sirona

Power Head  
24W (F/Optic)

HE21KL  
HE21NL  
HE21WL  
HE21SL

Small Head  
16W (F/Optic)

HE22KL  
HE22NL  
HE22WL  
- N/A -

Power Head  
24W (N/Optic)

HE21K  
HE21N  
HE21W  
- N/A -

Small Head  
16W (N/Optic)

HE22K  
HE22N  
HE22W  
- N/A -

**\$979 (FIBRE OPTIC)**

**\$693 (NON OPTIC)**

**BUY 3  
GET 10%  
OFF**

## ECO LINE CONTRA - ANGLES



Buy 3  
**GET 1 FREE**  
\*of equal or lesser value  
**Mix & Match**



- 1:5 Red
- Internal water
- For FG-burs
- (Ø 1.6 mm)

• LE15L (F/Optic)  
**\$1,595** **BUY 3 GET 1 FREE**

• LE15 (N/Optic)  
**\$1,375** **BUY 3 GET 1 FREE**

- 1:1 Blue
- Internal water
- For CA-burs
- (Ø 2.35 mm)

• LE11L (F/Optic)  
**\$1,155** **BUY 3 GET 1 FREE**

• LE11 (N/Optic)  
**\$891** **BUY 3 GET 1 FREE**

## BASIC LINE CONTRA - ANGLES



Buy 4  
Only  
**\$319 EA**  
**Mix & Match**



- 1:1 Blue
- Push button

• LB11/11  
**\$418** **BUY 3 \$339 EA**

**BUY 4 \$319 EA**

- 1:1 Blue
- For HP-burs

• LB02  
**\$418** **BUY 3 \$339 EA**

**BUY 4 \$319 EA**

## COUPLINGS



- Kavo type with LED and Water
- QC6016KW

**\$495**



- W&H type with Fibre Optic
- QC5016W

**\$385**



- NSK type with LED
- QC6016NT

**\$385**



- Sirona type with LED and Water
- QC6016SW

**\$495**

## AIR MOTORS

**NEW Shorter & Lighter**



- Internal water
- LED
- AM1006

**\$990**



- Internal water
- AM1004

**\$748**



- External water
- AM0014

**\$473**



- AM0004


**\$473**

## UNIVERSAL TORQUE WRENCH



**TORQ CONTROL®**  
For any implant

**3 Easy Steps**



Order code:  
15501

**\$1,210**

**BUY 2 GET 5% OFF**

**BUY 3 GET 10% OFF**

**Anthogyr**  
A Straumann Group Brand  
[www.anthogyr.com](http://www.anthogyr.com)

**Made in France**

**V ULTRASONIC SCALER**



**BUY 3  
GET 1  
FREE**



EMS Non-optic  
Type EN-041

- SC21E

\$462

**BUY 3  
GET 1 FREE**



EMS Piezon Non-optic  
Type EN-061

- SC21E2

\*Only for use with EMS compatible scaler tips

\$561

**BUY 3  
GET 1 FREE**



EMS Piezon LED  
Type EN-060A

- SC21EL

\$748

**BUY 3  
GET 1 FREE**



KaVo Piezo LED  
Type 1.007.3995

- SC21KL\*

\$748

**BUY 3  
GET 1 FREE**

**NEW**



Acteon-Satelec Non-optic  
Type F 12 281

- SC21S

\$462

**BUY 3  
GET 1 FREE**



Acteon-Satelec LED  
Type F 12 609

- SC21SL

\$748

**BUY 3  
GET 1 FREE**

**V PROPHY LINE AIR POLISHING DEVICE**

- Ergonomically shaped, compact and light-weight air-driven polishing device
- 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation



- PR1011

**\$1,375**

**Additional tip**



- SP3011

**\$329**

**Buy 1 Polisher**

Get 1 Bottle of Prophy Powder FREE

Get 50% DISCOUNT ON 1 ADDITIONAL TIP



- Prophy Powder (300 gms)  
Strawberry - PP1013  
\$ 60.50 each

Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N,  
W&H - PR1011W, Sirona - PR1011S

MANUFACTURER



**MK-DENT**

MK-dent GmbH  
Marie-Curie-Str. 2  
22941 Bargteheide  
GERMANY  
www.mk-dent.com

YOUR SPECIALIST DEALER

**ADAM DENTAL**  
supplies

Ph: 1300 449 262

• All catalogue prices are inclusive of GST • Promotion valid until 30th Apr 2025 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.